

❖ Tripp-Lite IV Pole Mount

- Securely mounts medical-grade power strip to IV pole; Track inventory & prevent loss
- Compliant with hospital standards
- Built-in splash guard helps protect power strip from liquid drips
- Cord hooks prevent damage and provide a safe, organized method for cord storage
- Eliminates trip hazards when not in use



Securely Mount Your Tripp-Lite Power Strip to IV Poles

The Tripp Lite IV Pole Mount was designed to securely mount your Tripp Lite Medical Grade Power Strip (415-HG OEM & PS-615-HG OEM) to any IV pole of up to two inches in diameter. Its robust cord wrap hooks aid with organization and will help to prevent cord damage or trip-hazards when not in use. Best of all, the Tripp Lite IV Pole Mount requires a an allen wrench tool in order to install/remove so power strip loss/theft is reduced.

Part Numbers: Manifold BlockPart Numbers: Available Flow Meters & DISS Fittings

Tripp Lite IV Pole Mount. Contact Alpha Souce